

Pacer

The Rifton Pacer's versatile design adapts to meet each person's unique needs, supporting users across all ability levels in any setting—from therapy clinics to home environments.



Scan here to order, request a quote, or see detailed product information.


- Promotes natural weight-bearing and weight-shifting
- Tool-free adjustments
- Choose between a standard base with advanced casters or utility base with large wheels.



Key features


1 Dynamic movement

The Rifton Pacer enables weight-bearing and weight-shifting movements necessary for a natural gait pattern.




2 Therapeutic positioning

The multi-position saddle (MPS) adjusts in five directions to accommodate each client's posture.




Video: Get to know the MPS adjustments.

The Pacer in use




The multi-position saddle (MPS) is the ultimate pelvic positioning aid.




The lightweight chest pad and handlebars configuration provides enhanced freedom of movement for users who need minimal support.

3 Standard or utility base

The standard base is designed for smooth surfaces. The utility base can be used outdoors.





Standard base
Casters: 5½"




Utility base
Front casters: 8"
Rear wheels: 11½"

4 Easy height adjustment

Adjust height in seconds—one button press unlocks the gas spring for smooth, effortless positioning.



The large-wheeled utility base glides smoothly over outdoor terrain.



The dynamic Pacer promotes weight-bearing and weight-shifting needed for a natural gait pattern.

Key features


5 Versatility

Pair an upper frame with a base frame. Then add the accessories that meet your client's specific needs.



6 Transportable


Press the button on the main column to detach the upper frame from the base for easy transport.



The Pacer in use



The odometer makes it easy to record client progress in any setting, in either feet or meters.




The attendant guide bar allows caregivers to guide and support users from in front or behind.

Pacer options




Scan for sizing and pricing.


Frame size options



Mini K610
Elbow Height:
15½ – 20½"




Small K620
Elbow Height:
18½ – 27½"



Medium K630
Elbow Height:
24 – 35"



Large K640
Elbow Height:
32 – 47"




XL K650
Elbow Height:
34 – 49"


Frame color options




Red



Blue




Black






Lime



Pink





Pacer components



Scan for sizing
and pricing.

Upper frame options




Dynamic upper frame: Allows weight-shifting and dynamic body-weight support, encouraging natural gait patterns.




Standard upper frame: A light-weight choice for clients who do not require dynamic movement.

Base options



Standard base


Designed to be used on smooth surfaces.



Utility base


This base is the best choice for navigating thresholds, lawns, gravel paths or chipped playgrounds.

Optional components




Arm prompts

Height, rotation, angle, and positioning adjustments.



Arm platforms

A simple alternative for users needing minimal arm positioning.



Hand brakes


Lever action provides instant stopping power and locks in place for parking.

More Pacer components



Hand loops

A simpler option for upper extremity support.




Chest prompt

The padded chest prompt holds the client's torso securely at any angle.



Chest pad


The single support strap and double angle adjustment allow for easy transfers and quick positioning.



Handlebars


Quick to attach, easy to adjust, and pair well with the chest pad.

Height from pivot to top: small: 5", large: 7"




Multi-position saddle

Adjusts in five directions for precise pelvic positioning.




Pelvic support

For users requiring minimal abduction. Assists weight-bearing while promoting freedom of movement and strength building.




Hip positioner

Enables natural pelvic movement and weight-shifting with good abduction. Optional comfort pad available.




Thigh prompts

Prevent leg scissoring and rotation within the frame. Adjustable positioning.




Ankle prompts

Prevent scissoring and control stride length with padded straps and spring adjusters.



Handholds


Ergonomic handholds adjust in height and position anywhere along the top bar. For use by caregiver or user.



Communication tray


Provides convenient access to tablets and mobile devices.

Inside dimensions: 8½" x 11"
Max. working load: 10 lb




Attendant guide bar

Attaches with a quick-release clamp, enabling caregivers to guide from front or behind.



Accessories tote

Keeps all Pacer components together. Large, durable, and attractive design.



Need support?

The Rifton website offers everything you need:

- Evidence-based articles
- Sample letters of medical necessity
- Video tutorials and training resources
- Product instruction guides
- Clinical resources and best practices



Scan here to order, request a quote, or see detailed product information.